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INSTRUCTIONS FOR USE BioFlote™ 4000 PNEUMATIC MATTRESS

NOTE!

Before using the mattress, please read the following instructions carefully.

I. INTRODUCTION

Dear Readers,

Thank you for choosing a product by REHA FUND Sp. z o.o., a leading supplier of medical supplies and rehabilitation products on the Polish market.

The Instructions for use contain tips and suggestions, the knowledge of which is essential for the correct device use. Keep this manual for future reference.

*Pictures in the instructions below are for illustrative purposes only.

GENERAL INFORMATION

The variable-pressure anti-bedsore mattress system uses a very advanced technology to ensure the highest quality device and comfort. It is designed for the prevention and treatment of bedsores.

The system has been tested and complies with the following standards:

EN 60601-1-11:2010, EN60601-1:2006, EN 60601-1-2: 2007, CISPR 11:2009, IEC 61000-4-3:2010, IEC 61000-4-4:2004, IEC 61000-4-5:2005, IEC 61000-4-6:2008, IEC 61000-4-8:2009, IEC 61000-4-11:2004, IEC 61000-3-2:2009, IEC 61000-3-3:2008.

EMC WARNING

This device has been tested and found compliant with limitations for medical devices EN 60601-1- 2:2007. These limitations are intended to provide adequate protection against harmful interference in a typical medical installation. This device generates, uses and may emit radio waves and, if not installed and used in accordance with the instructions, may cause interference to other devices in the vicinity. However, there is no guarantee that interference will not occur with a particular installation. If this equipment does cause interference to other equipment, which can be determined by turning the equipment on and off, the user may attempt to correct the interference by one or more of the following measures:

- change the orientation or position of the receiving device;
- increase the distance between the devices;
- connect the device to an outlet on a circuit different from that to which other devices are connected;
- ask the manufacturer, supplier or service technician for assistance.

NOTE!

The device is not suitable for use near flammable liquids and gases.

II. SPECIFICATIONS

TECHNICAL DATA

LIST	LIST OF PUMP SPECIFICATIONS
Model	BioFlote™ 4000 pneumatic mattress
Power supply	230 V/50 Hz
Fuse	1A
Pump dimensions in cm	28 x 15 x 9
Pump weight	1.4 kg
Power consumption	7 W
Cycle duration	12 min.
Pressure range	40-100 mmHg
Pump capacity	8 litre/min

LIST	LIST OF PUMP SPECIFICATIONS
Model	BioFlote™ 4000
Type	Pocket (17 air pockets)
Dimensions of the mattress in cm	190 x 88 x 11 – after inflation
Weight of the mattress	4.85 kg
Maximum user weight	140 kg
Material	Nylon/Polyvinyl chloride
Colour	blue
Ambient temperature	Working: 5°C to 40°C Storage: -25°C to 70°C Transport: -25°C to 70°C
Ambient relative humidity:	Working: 15%~90% Storage: 10%~90% Transport: 10%~90%
Atmospheric pressure	40–100 mmHg

NOTE

1. The above technical parameters are valid at all sites with the same electricity parameters.
2. For other technical documents, contact your supplier.
3. Dispose the product in accordance with local legal requirements.

APPENDIX A: EMC INFORMATION

Guidelines and manufacturer's declaration – electromagnetic emissions

This device is intended for use in the electromagnetic environment described below. The user of this device should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidelines
RF emissions CISPR 11	Group 1	The device uses RF energy exclusively for internal functions. Therefore, RF emissions are very low and cannot cause any interference to nearby electronic devices.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000–3–2	Class A	The device is suitable for use in all rooms, including domestic premises, and can be connected to the public low-voltage power grid.
Voltage fluctuations/ vibration emissions IEC61000–3–3	Complies	

Guidelines and manufacturer's declaration – Electromagnetic Resistance

This device is intended for use in the electromagnetic environment described below. The user of this device should ensure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) IEC61000–4–2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC61000–4–4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	The quality of the power supply should be typical (not as for industry or hospitals)

Wave IEC61000–4–5	± 1 kVlines ± 2 kVlines – ground	± 1 kVlines	The quality of the power supply should be typical (not as for industry or hospitals)
Voltage dips, short interruptions and voltage fluctuations on power supply input lines	<5%UT (>95% spike in UT) for 0.5 cycles 40%UT (60% spike in UT) for 5 cycles 70%UT (30% spike in UT) for 25 cycles <5%UT (>95% spike in UT) for 5 seconds		The quality of the power supply should be typical (not as for industry or hospitals) If the user wants uninterrupted and undisturbed operation of the device, it should be provided with some permanent power supply or batteries.
Power frequency (50/60Hz) magnetic field IEC6100–4–8	3 A/m	3 A/m	The frequency strength of the magnetic field should be at the typical level for atypical locations, in a typical industrial or hospital environment
NOTE: UT is the value of the AC supply voltage before testing			

Guidelines and manufacturer’s declaration – Electromagnetic Resistance

This device is intended for use in the electromagnetic environment described below. The user of this device should ensure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance	Electromagnetic environment – guidelines
RF conductivity IEC61000–4–6	3 Vrms 150 kHz to 80 MHz outside the ISM band ^a	3Vrms	Portable and mobile RF devices should be used no closer to any part of the device, including cables, than the recommended distance calculated from the equation appropriate to the frequency of the transmitter. Recommended distance: $d = [3 \frac{\sqrt{P}}{4\pi f}]$ $\sqrt{P}$ from 150 kHz to 80 MHz $d = [3 \frac{\sqrt{P}}{4\pi f}]$ $\sqrt{P}$ from 80 MHz to 800 MHz $d = [7 \frac{\sqrt{P}}{4\pi f}]$ $\sqrt{P}$ from 800 MHz to 2.5 GHz where P is the maximum power rating of the transmitter in watts (W) according to the transmitter manufacturer, and d is the recommended distance in metres (m). ^b The field strength from fixed RF transmitters, as determined by the electromagnetic site survey ^c , should be less than the compliance level at each fluctuation frequency. d Interference may occur in the vicinity of equipment marked with the following symbol: 
RF radiation IEC61000–4–3	3 V/m from 80 MHz to 2.5 GHz	3 V/m	


PNEUMATIC PUMP CONSTRUCTION

- Quick connector;
- Front control panel;
- Transport connector;
- Power cord;
- Air filter;
- Pump attachments.


FRONT PANEL

1. Air pressure regulator – is responsible for the air pressure value in the mattress. Turning the knob towards "MAX" increases the pressure, turning the knob towards "MIN" decreases the pressure in the mattress.
2. Main power switch – to switch the device's power on and off.

CONSTRUCTION OF THE ANTI-BEDSORE KIT

1. CPR valve;
2. Pressure-controlled pump;
3. Pocket mattress with cover;
4. Quick connector;
5. Transport connector.

III. PURPOSE AND USE OF THE PRODUCT

This product is designed for people:

- in need of treatment and prevention of grade II bedsores while optimising patient comfort;
- permanently in a home or hospital position, long-term care and treatment of bedsores;
- in need of treatment and reduction of pain caused by ground pressure on the body's soft tissue when prescribed by a doctor or physiotherapist.

This product may be operated by people who are qualified to perform general nursing procedures and have the appropriate training and knowledge of pressure ulcer prevention and treatment.

IV. CONTRAINDICATIONS

Do not use in patients with spinal cord injury.

V. GENERAL RULES FOR SAFE USE

SAFETY RULES

When using electrical products, and especially in the presence of children, all safety measures, including the following, must be observed.

DANGER!

To reduce the risk of electric shock:

1. Do not use the device while bathing.
2. Always switch off the device immediately after use.
3. Do not place or store the device where it can fall from or be pulled into a bathtub or washbasin.
4. Do not place or immerse the device in water or any other liquid.
5. Do not touch a device that has fallen into water. If this happens, immediately unplug the device.

WARNINGS!

To reduce the risk of electric shock, fire or injury:

1. Assess the likelihood of the device becoming a patient trap and limiting the patient's range of movement or trapping them completely. If this is likely to occur, the patient should be kept under continuous monitoring.
2. Do not use in patients with spinal cord injury.
3. It is essential to supervise the performance of the device when it is used by or near children. Some small components can be easily dismantled, posing a risk of choking or electrocution to children.
4. Use for the purpose described in the instructions for use. Do not use for purposes not recommended by the manufacturer. Do not use any type of mattress other than that indicated by the manufacturer.
5. Do not use if the device is not working properly, the cord or plug is damaged, has been dropped, damaged or has fallen into water. In this case, the device should be taken to an authorised service centre for inspection/repair.
6. Keep the power cord away from hot surfaces.
7. Do not obstruct the air intake. Do not place on a soft surface (bed, mattress, blanket) which may block the air inlet. Keep the immediate area around the pump clean – away from hair, dander, clippings or other small particles.


8. Do not modify/repair the device yourself without the manufacturer's authorisation.
9. If you suspect that you may have or have already had allergic symptoms after contact with the mattress, contact your doctor immediately.
10. Each serious incident related to the product shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient reside.

VI. METHOD OF USE

INSTALLATION AND UNPACKING

The pump and mattress are packed in a cardboard box to protect the contents inside. Remove the packaging and check that nothing has been damaged in transit. In the event of damage, contact your dealer immediately.

SETUP

1) Arrange the mattress on the bed. It is recommended to arrange the mattress so that the end of the mattress with the inlet ends is on the feet side. Ensure that the flaps on both sides of the mattress are tucked under the bed mattress.

NOTE!

The mattress should be covered with a sheet to prevent skin contact with the mattress in order to improve patient comfort, and to protect the mattress from contact with moisture or medicines (creams, ointments).

2) Unfold the hooks and hang the pump on the edge of the bed or place it on the floor in a suitable location to avoid exposure of the pump to damage.

3) Connect the air tubes to the pump and mattress. A typical "click" will indicate that the wires are correctly connected.

NOTE!

Check and ensure that the air tubes are not twisted or pinched.

4) Connect the pump with the power cord to the socket.

NOTE!

Before inserting the plug into the socket, ensure that the power supply is compatible with the device's requirements. Using the power plug, the device can also be disconnected. The electrical socket should be in such a position that it allows the device to be quickly unplugged.

5) Switch on the pump power supply (switch to ON position).

6) Wait approximately 20 minutes for one of the sections to inflate completely before placing the patient on the mattress. Set the regulator knob to maximum and, once the mattress is inflated, set the pressure to suit the patient's weight.

NOTE!

The pump may only be used with the manufacturer's recommended mattress. Do not use the pump for other purposes. It is important to note that the pump must be running at all times when the mattress is in use. Once the pump has been switched off, the mattress should not be used. Do not pump the mattress while the patient is lying on it.


A few comments on the installation:

Once set up, make sure that the power cable lies so that it does not pave the way or cause tripping or falling. The device should be positioned so that the doctor or user has easy access to the device and can disconnect it at any time.

FUNCTIONS AVAILABLE

◆ CPR valve

If for any reason there is a need to quickly deflate the mattress, pull the CPR valve handle. The CPR valve is located at the head end of the mattress.


◆ Pressure adjustment

Using the pressure regulator, the user can set the desired softness of the mattress. The setting of the firmness level of the mattress should be consulted with your GP.

◆ Transport

A quick-release coupling can be used to allow the inflated mattress to be moved to block the air outflow from the mattress for transport.

TROUBLESHOOTING

NOTE!

To prevent the risk of electric shock, do not attempt to remove the housing. Any repairs must be reported to an authorised service centre. The included repair kit can be used after the warranty period has expired. Repairing the mattress without reporting the damage to an authorised service centre will result in immediate loss of warranty rights.

PROBLEM	DISSOLUTION
No power	Check that the plug is inserted into an outlet where there is voltage.
Low pressure indicator light is lit	Check that the air lines to the pump and mattress are correctly connected. Check the connections between the pockets. Check for leaks at the CPR valve. If the pump has a pressure control, increase the pressure. Check that the mattress pockets are not perforated or cracked. NOTE: The pump must be switched on at all times when using the mattress.
A noise can be heard from time to time in the pump	This is a normal symptom. When one section is emptied, air escapes through the pump, which can be heard.
Patient collapses	The pressure value set on the pump with the regulator is not suitable for the patient. Increase the pressure and wait for the air in the pockets to refill.
The mattress does not hold its shape	Check that the pockets are held in place with straps attached to the base of the mattress, that the pockets are fastened with snaps to the base.
Part of the mattress is uninflated	This is a normal phenomenon when one section is inflated and the other is not. After half the cycle has elapsed, the pumping of the section should change. If only one of the sections is pumping at all times, this may indicate a leaking mattress or damage to the pump.

If the above advice has not helped to resolve the fault, please contact your local representative directly. You will probably need to consult a technical service.

VII. CLEANING AND MAINTENANCE

PUMP

Clean the pump with a damp cloth and mild detergent and keep away from dust. Detergents that will cause a chemical reaction on the surface of the plastic housing of the pump must not be used. All parts should be thoroughly dried before re-use.

NOTE!

Ensure that the pump is disconnected from the power supply before cleaning. Do not immerse or wet the pump.

MATTRESS

Clean the mattress by wiping with a cloth or sponge moistened in warm water with a mild detergent. All parts should be thoroughly dried before re-use.

NOTE!

Do not use solvent- or alcohol-based detergents to clean the mattress. This may cause damage to the surface of the mattress in the form of cracks or hardening of the material. After cleaning, dry the mattress in a shaded area.

STORAGE

1. Place the deflated mattress on a flat surface.
2. Roll the mattress from head to foot towards the air outlet ends with the CPR valve open.
3. Secure the pipes by inserting them into the mattress.

NOTE!

The mattress must not be folded or crushed. Avoid direct contact with the sun.

INSPECTION

GENERAL

1. Check for wear on the power cord and plug.
2. Check the surface of the mattress for wear, hardness or damage. Ensure that the pockets of the mattress are correctly connected to each other, check the tightness of the pockets.
3. Check for damage to the duct connections and to the air ducts themselves (kinks, cracks).
4. Check the correct performance of the pump. The airflow should change between the outlets every 1/2 cycle period.

AIR FILTER REPLACEMENT

The air filter should be checked and replaced if necessary (every 12 months or so). The fouling of the air filter depends on the environment in which the pump is operating. A dirty (clogged) filter will cause the compressor inside the pump to malfunction and even damage it.

1. Open the air filter cover located at the rear of the pump housing by releasing the catch.
2. Check air filter and replace if necessary.
3. Fit the filter cover.

VIII. END OF USE

If abnormalities are found during the mattress inspection in accordance with the manufacturer's recommendations. See section: Cleaning and Maintenance – stop using the product immediately. At the end of use, the product must be disposed of in accordance with local environmental regulations.

IX. INTENDED PERIOD OF USE

The product ensures safe and reliable performance during the warranty period if the product is used and maintained in accordance with the manufacturer's instructions for use'. The manufacturer recommends that the kit is only checked and serviced by an authorised service provider if there are signs of wear or suspected faulty operation of the device. In other cases, service or inspection of the device is generally not required.

NOTE!
A spare pocket is included in the kit. If the mattress is found to be leaking, contact an authorised service centre.

X. WARRANTY:

The manufacturer assures the customer that the product sold is brand new and of high quality. The product is covered by a warranty (see the manufacturer's website for detailed warranty conditions). The product should be used as intended and as recommended in the Instructions for use.

XI. GLOSSARY OF SYMBOLS

	Manufacturer		Top, do not roll over
	CE marking		Protect from humidity
	Catalogue No.		Fragile, handle with care
	Date of manufacture		Type BF parts
	See the instructions for use		Protection class I device
	Warning	IP22	Protection against solid objects 12.5mm and larger. Protection against vertically falling water drops
	Batch code		The device is marked in accordance with the European Directive (WEEE) and the Polish Act on Waste Electrical and Electronic Equipment. This marking indicates that the equipment must be returned to a waste electrical and electronic equipment treatment facility at the end of its useful life.
	Medical device		Do not cut
	Unique product identification code		Recycling

REHA FUND Sp. z o.o.
ul. Staniewicka 14,
03–310 Warsaw, Poland
www.rehafund.pl

I/BioFlote4000/3/08.11.2024
Instructions for use issue date: 11.2024